

Osteopath Laser Phototherapy System

Xp Control Unit

Choose between either the softly moulded Clinic unit in white with matching probes or the compact grey Mobile unit with carry bag. Both have built-in beam test, energy density calculator, long-life rechargeable battery for non-mains use and a trigger point locator function with outcome change measurement.

Infra-Red Laser Cluster

For subcutaneous soft tissue repair and pain relief across larger treatment areas we recommend the 5 x 820nm 200mW infra-red cluster probe. The 820nm wavelength is a real all-rounder with good response across a range of relevant cells and capable of tackling chronic adhered tissue as well as optimising resolution of acute musculoskeletal injuries.

Large Multi-Wavelength LED Cluster

This 60 diode cluster is excellent for reducing oedema and its set of 5 wavelengths work synergistically in combination with the single wavelength probes, being particularly helpful for consolidating adjustment and improving client comfort post-manipulation.

Infra-Red Single Laser

Indispensable for trigger point application and non-needle acupuncture we recommend inclusion of the single diode 820nm 200mW. This is also effective for closely targeted tissue repair and pain relief.

Point-Finder Kit

Each Xp Control Unit comes with a dedicated 'hand probe' and a specialized tip for a single laser probe to enable the point-locator function.

On-site training

Product Introductory training at the customer's site for up to 6 people is provided within the package price (for mainland UK direct purchases).

Two Pairs of Protective Glasses

Eye protection appropriate to the wavelengths utilised for use by clinician and patient.



www.omegalaser.co.uk

Osteopath System
with Xp Mobile Control Unit

Professional Quality

Omega equipment is used on a truly worldwide basis, with systems used in over 70 countries. Our latest products have had USA FDA clearance, are fully licensed for Canada and have all relevant European certification for medical use.

Training

We recognise the benefits to all in ensuring that the laser system fits the clinician's needs and that they can employ it effectively, so we advise on appropriate tailored systems and train one-to-one with each direct purchase.

After-Sales Service

Our strong emphasis on after-sales service distinguishes us in the field with efficient and responsive processes plus ongoing phone support for new customers as they build confidence in using their equipment for the host of treatments possible.

The Company

Omega Laser Systems Limited is the longest established British manufacturer of medical standard Low Level Laser Equipment. It is a specialist in 3B lasers with some 25 years of research, development and manufacturing experience in this field. The laser systems are made and serviced at Omega's UK workshops.

Flexible Systems

An Omega Laser System comprises a control unit plus a selection of interchangeable probes according to the treatments to be undertaken. This leaflet describes the system recommended for Osteopaths but do get in touch if you would like to discuss whether an alternative system would better suit your needs or budget. For more information on equipment refer to our general brochure, our website or call us on 01787 477551.



Osteopath System
with Xp Clinic Control Unit

www.omegalaser.co.uk